



ADULT STETHOSCOPE (SD5030 Series)



SPECIFICATIONS:

- Traditional bell
- 47mm diaphragm
- Product's versatility provides clinicians added options & increases level of care
- Bell side = low frequency (lung)
- Diaphragm = high frequency (heart)
- 1 stethoscope per box
- 100 stethoscopes per case
- Box and Case are UPC bar coded

FEATURES/BENEFITS:

- Ergonomic brass binaural
- Non Latex
- Lightweight metal chest piece with a non-chill rim
- 24" vinyl tubing
- Single patient use design helps reduce the risk of cross-contamination
- Several different color options can be matched to a facility's standardized color system

